

## Reaction of Sodium Thiosulphate with Hydrochloric acid

**Safety:** Good laboratory practice should be followed. Concentrations of hydrochloric acid stronger than 2M but weaker than 6.5M should be labelled IRRITANT but this investigation does not need concentrations approaching 2M

**Reaction:** Chloride displaces thiosulphate. In the first reaction thiosulphuric acid is formed. This is unstable and decomposes producing sulphur dioxide and solid sulphur that is forms the cloudy precipitate.





1.  $\text{Na}_2\text{S}_2\text{O}_3 (\text{aq}) + 2\text{HCl} (\text{aq}) \rightarrow 2\text{NaCl} (\text{aq}) + \text{H}_2\text{S}_2\text{O}_3(\text{aq})$
2.  $\text{H}_2\text{S}_2\text{O}_3(\text{aq}) \rightarrow \text{S} (\text{s}) + \text{SO}_2 (\text{g}) + \text{H}_2\text{O} (\text{l})$

The rate of the reaction depends on the concentrations of acid and sodium thiosulphate and on the temperature. These dependencies can be demonstrated and accurately measured using the colorimeter.

### Using the Mystrica colorimeter to measure reaction rates:

Measure the absorbance of blue light. The colorimeter can be used stand-alone or connected to a computer.





See a video clip of this reaction on the *Mystrica* website

[www.mystrica.com](http://www.mystrica.com)