

0.1M Citric acid-Sodium citrate buffer buffer – pH range 3.0 – 6.2

Prepare a 0.1M solution of citric acid monohydrate, C₆H₈O₇H₂O (21.01g/l) and a 0.1M solution of trisodium citrate dihydrate, C₆H₅O₇Na₃·2H₂O (29.41g/l). Mix the volumes shown in the table.

Or dissolve the masses shown and make up to 100cm³ with water

pH	0.1M citric acid (cm ³)		0.1M trisodium citrate (cm ³)	
	Volume of 0.1M solution (cm ³)	Mass in 100cm ³ (g)	Volume of 0.1M solution (cm ³)	Mass in 100cm ³ (g)
3.0	82	1.72	18	0.53
3.2	77.5	1.63	22.5	0.66
3.4	73	1.53	27	0.79
3.6	68.5	1.44	31.5	0.93
3.8	63.5	1.33	36.5	1.07
4.0	59	1.24	41	1.21
4.2	54	1.13	46	1.35
4.4	49.5	1.04	50.5	1.49
4.6	44.5	0.93	55.5	1.63
4.8	40	0.84	60	1.76
5.0	35	0.74	65	1.91
5.2	30.5	0.64	69.5	2.04
5.4	25.5	0.54	74.5	2.19
5.6	21	0.44	79	2.32
5.8	16	0.34	84	2.47
6.0	11.5	0.24	88.5	2.60
6.2	8	0.17	92	2.71